

# Use of force-controlled grippers

## 1 Python control API description

**Note** : The time interval between calling each instruction of the gripper must be greater than 1.5 seconds.

### `set_pro_gripper(gripper_id, address, value)`

- **Function** : Set the parameters of the Pro force-controlled gripper. You can set a variety of parameter functions. Please see the following table for details.
- parameter:
  - `gripper_id` ( `int` ): gripper ID, default value is 14, value range is 1 ~ 254.
  - `address` ( `int` ): The command number of the gripper.
  - `value` : the parameter value corresponding to the instruction number.

| Function                                             | <code>gripper_id</code> | <code>address</code> | <code>value</code>              |
|------------------------------------------------------|-------------------------|----------------------|---------------------------------|
| Set the gripper ID                                   | 14                      | 3                    | 1 ~ 254                         |
| Set the gripper enable state                         | 14                      | 10                   | 0 or 1, 0 - disable; 1 - enable |
| Set the clockwise runnable error of the gripper      | 14                      | 21                   | 0 ~ 16                          |
| Set the anti-clockwise runnable error of the gripper | 14                      | 23                   | 0 ~ 16                          |
| Setting the minimum actuation force of the gripper   | 14                      | 25                   | 0 ~ 254                         |
| IO output settings                                   | 14                      | 29                   | 0, 1, 16, 17                    |
| Set the IO opening angle                             | 14                      | 30                   | 0 ~ 100                         |
| Set IO closing angle                                 | 14                      | 31                   | 0 ~ 100                         |
| Set the virtual position value of the servo          | 14                      | 41                   | 0 ~ 100                         |
| Setting the holding current                          | 14                      | 43                   | 1 ~ 254                         |

- Return Value
  - Please see the following table:

| Function                     | Return Value               |
|------------------------------|----------------------------|
| Set the gripper ID           | 0 - failed; 1 - successful |
| Set the gripper enable state | 0 - failed; 1 - successful |

| Function                                             | Return Value               |
|------------------------------------------------------|----------------------------|
| Set the clockwise runnable error of the gripper      | 0 - failed; 1 - successful |
| Set the anti-clockwise runnable error of the gripper | 0 - failed; 1 - successful |
| Setting the minimum actuation force of the gripper   | 0 - failed; 1 - successful |
| IO output settings                                   | 0 - failed; 1 - successful |
| Set the IO opening angle                             | 0 - failed; 1 - successful |
| Set IO closing angle                                 | 0 - failed; 1 - successful |
| Set the virtual position value of the servo          | 0 - failed; 1 - successful |
| Setting the holding current                          | 0 - failed; 1 - successful |

### get\_pro\_gripper(gripper\_id, address)

- **Function** : Get the parameters of Pro force-controlled gripper, and you can get various parameter functions. Please see the following table for details.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
  - `address ( int )`: The command number of the gripper.

| Function                                              | gripper_id | address |
|-------------------------------------------------------|------------|---------|
| Read the firmware major version number                | 14         | 1       |
| Read the firmware minor version number                | 14         | 2       |
| Read the gripper ID                                   | 14         | 3       |
| Read the clockwise runnable error of the gripper      | 14         | 22      |
| Read the anti-clockwise runnable error of the gripper | 14         | 24      |
| Read the minimum actuation force of the gripper jaws  | 14         | 26      |
| Read IO opening angle                                 | 14         | 34      |
| Read IO closing angle                                 | 14         | 35      |
| Get the amount of data in the current queue           | 14         | 40      |
| Read the virtual position value of the servo          | 14         | 42      |
| Reading the holding current                           | 14         | 44      |

- Return value:
  - Check the following table (if the return value is -1, it means no data can be read):

| Function                                              | Return Value                                                     |
|-------------------------------------------------------|------------------------------------------------------------------|
| Read the firmware major version number                | Major version number                                             |
| Read the firmware minor version number                | Minor version number                                             |
| Read the gripper ID                                   | 1 ~ 254                                                          |
| Read the clockwise runnable error of the gripper      | 0 ~ 254                                                          |
| Read the anti-clockwise runnable error of the gripper | 0 ~ 254                                                          |
| Read the minimum actuation force of the gripper jaws  | 0 ~ 254                                                          |
| Read IO opening angle                                 | 0 ~ 100                                                          |
| Read IO closing angle                                 | 0 ~ 100                                                          |
| Get the amount of data in the current queue           | Returns the amount of data in the current absolute control queue |
| Read the virtual position value of the servo          | 0 ~ 100                                                          |
| Reading the holding current                           | 1 ~ 254                                                          |

### set\_pro\_gripper\_angle(gripper\_id, gripper\_angle)

- **Function** : Set the angle of the force-controlled gripper.
- parameter:
  - `gripper_id` ( `int` ): gripper ID, default value is 14, value range is 1 ~ 254.
  - `gripper_angle` ( `int` ): Gripper angle, value range 0 ~ 100.
- Return value:
  - 0 - Failed
  - 1 - Success

### get\_pro\_gripper\_angle(gripper\_id)

- **Function** : Read the angle of the force-controlled gripper jaws.
- parameter:
  - `gripper_id` ( `int` ): gripper ID, default value is 14, value range is 1 ~ 254.
- **Return value** : `int` 0 ~ 100

### set\_pro\_gripper\_open(gripper\_id)

- **Function** : Open the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:
  - 0 - Failed
  - 1 - Success

### set\_pro\_gripper\_close(gripper\_id)

- **Function** : Close the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:
  - 0 - Failed
  - 1 - Success

### set\_pro\_gripper\_calibration(gripper\_id)

- **Function** : Set the zero position of the force-controlled gripper. (The zero position needs to be set first when using it for the first time)
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:
  - 0 - Failed
  - 1 - Success

### get\_pro\_gripper\_status(gripper\_id)

- **Function** : Read the clamping status of the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:
  - 0 - Moving.
  - 1 - Motion stopped, no clamping detected.
  - 2 - Stop motion, clamping object detected.
  - 3 - After the clamping object is detected, the object falls.

### set\_pro\_gripper\_torque(gripper\_id, torque\_value)

- **Function** : Set the torque of the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.

- `torque_value ( int )`: Torque value, ranging from 1 to 100.
- Return value:
  - 0 - Failed
  - 1 - Success

### `get_pro_gripper_torque(gripper_id)`

- **Function** : Read the torque of the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
- **Return value:** ( `int` ) 100 ~ 300

### `set_pro_gripper_speed(gripper_id, speed)`

- **Function** : Set the force-controlled gripper speed.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
  - `speed (int)`: Gripper movement speed, ranging from 1 to 100.
- Return value:
  - 0 - Failed
  - 1 - Success

### `get_pro_gripper_default_speed(gripper_id, speed)`

- **Function** : Read the default speed of the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
- **Return value** : the default movement speed of the gripper, ranging from 1 to 100.

### `set_pro_gripper_abs_angle(gripper_id, gripper_angle)`

- **Function** : Set the absolute angle of the force-controlled gripper.
- parameter:
  - `gripper_id ( int )`: gripper ID, default value is 14, value range is 1 ~ 254.
  - `gripper_angle ( int )`: Gripper angle, value range 0 ~ 100.
- Return value:
  - 0 - Failed
  - 1 - Success

### `set_pro_gripper_pause(gripper_id)`

- **Function** : Pause movement.
- parameter:
  - `gripper_id ( int )`: Gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:

- 0 - Failed
- 1 - Success

### **set\_pro\_gripper\_resume(gripper\_id)**

- **Function** : Sports recovery.
- parameter:
  - `gripper_id ( int )` Gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:
  - 0 - Failed
  - 1 - Success

### **set\_pro\_gripper\_stop(gripper\_id)**

- **Function** : Stop motion.
- parameter:
  - `gripper_id ( int )` Gripper ID, default value is 14, value range is 1 ~ 254.
- Return value:
  - 0 - Failed
  - 1 - Success

## **2 Python zero calibration instructions**

After attaching the gripper to the end of the robotic arm, run the following code, and the gripper will first open to its maximum, then fully close to complete the calibration.

```
from pymycobot import MyCobot320
mc=MyCobot320("COM1")#Enter the serial port number of the robot
mc.set_pro_gripper_calibration(14)#The default ID of the gripper is 14
```

## **3 Python control scenario examples**

### **demo1: torque value adaptation**

The user can observe the deformation degree of the force-controlled gripper after clamping the object under the premise of setting the corresponding torque value and select a torque value with the smallest deformation but can clamp the object normally. This torque value is the most suitable torque value for the current object.

```

from pymycobot import MyCobot320
import time
mc=MyCobot320("COM1")#Fill in the serial port number of the robot
mc.set_pro_gripper_torque(14,10)#Set the torque value
print(mc.get_pro_gripper_torque(14))#View the set torque value
mc.set_pro_gripper_open(14)#Gripper fully opens
time.sleep(3)
mc.set_pro_gripper_close(14)#Gripper fully closes

```

## demo2: grab and move objects at points AB

```

from pymycobot import MyCobot320
import time
mc=MyCobot320("COM1")#Fill in the serial port number of the robot
mc.set_pro_gripper_torque(14,10)#Set the torque value of the gripper
print(mc.get_pro_gripper_torque(14))#Print the set gripper torque value
mc.set_pro_gripper_open(14)#First open the gripper completely
start_angles=[19.68, -1.23, -91.4, -0.52, 90.08, 60.29]
target_coords=[[231.3, -61.3, 232.7, 178.35, -2.7, -130.56],[231.3, 65.3, 232.7, 178.35, -2.7, -130.56]]
end_angles=[80,0,-85,0,90,60]
for i in range(len(target_coords)):
    mc.sync_send_angles(start_angles,50)
    mc.send_coords(target_coords[i],100,1)
    time.sleep(2)
    mc.send_coord(3,165,50)
    time.sleep(2)
    mc.set_pro_gripper_close(14)
    time.sleep(2)
    mc.send_coords(target_coords[i],100,1)
    mc.sync_send_angles(end_angles,100)
    mc.set_pro_gripper_open(14)
    time.sleep(2)

```

